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
**Table 2.** Average values of corneal riboflavin concentration ( $\pm$ SD) during the dosing phase of Thera-CXL, mean corneal stiffness parameter ( $k_c$ ,  $\pm$ SD), maximum deformation amplitude (DA,  $\pm$ SD) before and after Thera-CXL protocols, and average theranostic score in both protocols ( $\pm$ SD)

Parameter	Time (min)	Thera-CXL epi-on protocol (n = 12)	Thera-CXL epi-off protocol (n = 8)
Corneal riboflavin concentration ( $\mu\text{g}/\text{cm}^3$ )	5	50 $\pm$ 19	344 $\pm$ 146
	10	73 $\pm$ 36	558 $\pm$ 235
	15	118 $\pm$ 68	—
	20	145 $\pm$ 77	—
Mean corneal stiffness ( $k_c$ , N/m)	Baseline $k_c$	41.4 $\pm$ 12.7	34.0 $\pm$ 5.0
	After treatment $k_c$	49.9 $\pm$ 11.2 ( $P = .03$ )*	41.7 $\pm$ 7.0 ( $P = .04$ )*
Maximum deformation amplitude (DA, $\mu\text{m}$ )	Baseline DA max	635 $\pm$ 118	673 $\pm$ 86
	After treatment DA max	561 $\pm$ 104 ( $P = .03$ )*	610 $\pm$ 76 ( $P = .04$ )*
Theranostic score (d.n.)		0.7 $\pm$ 0.3	1.0 $\pm$ 0.2

d.n. = dimensionless number; Thera-CXL = theranostic-guided corneal crosslinking  
\*Statistically significant (paired t test)

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# STUDIO ARGO: LA CONVALIDA

ASSESSMENT OF THERANOSTIC GUIDED RIBOFLAVIN/UV-A CORNEAL

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Predictive ability of the theranostic imaging biomarkers in assessing the propensity of CXL treatment in flattening the  $K_{max}$  value at 1 year follow-up visit.

Predicted outcome	$K_{max}$ change at 1-year >0.1 D flattening	$K_{max}$ change at 1-year ( $\leq -0.1$ D flattening)	Total
Theranostic biomarkers (cut-off values)			
Riboflavin score >0.40 & Theranostic score $\geq 0.60$	True Positives (n = 41) $K_{max}$ change = $-1.7 \pm 1.6$ D	False Positives (n = 2) $K_{max}$ change = $+0.1 \pm 0.01$ D (epi-off only)	43
Riboflavin score $\leq 0.40$ and/or Theranostic score <0.60	False Negatives (n = 2) $K_{max}$ change = $-0.5 \pm 0.2$ D (epi-on only)	True Negatives (n = 2) $K_{max}$ change = $+0.2 \pm 0.2$ D (epi-on only)	4
<b>Total</b>	<b>43</b>	<b>4</b>	<b>47</b>

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
CROSS-LINKING FOR TREATMENT OF KERATOCONUS (REG.N.NCT05457647)<sup>19,20</sup>

## $K_{MAX}$ CHANGE IN TRUE POSITIVES OF STRATIFICATION GROUPS (1 year)

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
Regensight is a Vision Engineering Italy's Company

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## Dicono di noi

**Excellence in sight: trusted evidence, informed decisions, and better vision health**

Join us in ARVO Special Interest Group (SIG) Virtual Meeting Room  
Thursday, May 12, 2022  
2:00-3:30 PM MDT

Trusted evidence. Informed decisions. Better health.

### CONCLUSIONS of META-ANALYSIS STUDY

"Most effective method for **corneal cross-linking** is not known. **Theranostic UV-A device** measuring real-time corneal concentration of riboflavin, via acquisition and analysis of fluorescence emitted from riboflavin when illuminated by UV-A light (*Marco Lombardo et al.*) **to solve it**".

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EU Startup News (Novembre 2023)

## Which Italian Pharmaceutical Startups Are Shaping Europe's MedTech Future?

by EU Startup News

### Regensight

Regensight focuses on enhancing human vision through incision-free theranostics technology. A blend of therapeutics and diagnostics, their innovative solutions represent a big leap in vision care.

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### Bibliografia

1. Lombardo M et al. Translational label-free nonlinear imaging biomarkers to classify the human corneal microstructure. *Biomed Opt Express* 2015; 8;6(8):2803-2818
2. Alizadeh M et al. Identifying crossing collagen fibers in human corneal tissues using pSHG images. *Biomed Opt Expr* 2019;10(8):3875-3888
3. Lombardo M et al. Estimated prevalence of keratoconus in the largest metropolitan area of Italy *Eur J Ophthalmol* 2024; 34(3):649-655
4. Orphanet. Isolated keratoconus non rare disease (<https://www.orpha.net/>)
5. Ferdi AC et al. Keratoconus Natural Progression: A Systematic Review and Meta-analysis of 11 529 Eyes. *Ophthalmology* 2019; 126(7):935-945
6. Antunes-Foschini A et al. Natural history and predictors for progression in pediatric keratoconus. *Sci Rep* 2023; 13(1): 4940
7. Sykasis E et al. Corneal collagen cross-linking for treating keratoconus. *Cochrane Database Syst Rev.* 2015;2015(3):CD010621
8. Ng SM et al. Transepithelial versus epithelium-off corneal crosslinking for progressive keratoconus. *Cochrane Database Syst Rev* 2021;3(3):CD013512
9. Serrao S et al. Adverse events after riboflavin/UV-A corneal cross-Linking: a literature review. *Int Ophthalmol* 2022; 42(1): 337-348
10. Lombardo M et al. Real time monitoring of riboflavin concentration using different clinically available ophthalmic formulations for epi-off and epi-on corneal cross-linking. *Graefes Arch Clin Exp Ophthalmol* 2024; 262: 2569-2577
11. Chungyu T et al. Corneal collagen cross-linking in keratoconus: a systematic review and meta-analysis. *Sci Report* 2014; 4: 5652
12. Li J et al. Efficacy of corneal collagen cross-linking for treatment of keratoconus: a meta-analysis of randomized controlled trials. *Plos ONE* 2015;10(5):e0127079
13. Lombardo G et al. Non-invasive optical method for real-time assessment of intracorneal riboflavin concentration and efficacy of corneal cross-linking. *J Biophotonics* 2018;11(7):e201800028

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### Bibliografia

14. Lombardo M et al. Non-invasive and real time assessment of riboflavin consumption in standard and accelerated corneal cross-linking. *J Cataract Refract Surg* 2019;45(1):80-86
15. Lombardo G et al. Comparison between standard and transepithelial corneal cross-linking using a theranostic UV-A device. *Graefes Arch Clin Exp Ophthalmol* 2020;258(4):829-834
16. Lombardo G et al. Theranostic-guided corneal cross-linking: pre-clinical evidence on a new treatment paradigm for keratoconus. *J Biophotonics* 2022;15(12):e202200218
17. Lombardo M et al. Predicting corneal cross-linking treatment efficacy with real-time assessment of corneal riboflavin concentration. *J Cataract Refract Surg* 2023; 49(6): 635-641
18. Lombardo M et al. Spatial targeted delivery of riboflavin with a controlled corneal iontophoresis delivery system in theranostic-guided UV-A light phototherapy. *J Biophotonics* 2024; 17(7):e202400068
19. Roszkowska AM et al. A randomized clinical trial assessing theranostic-guided corneal cross-linking for treating keratoconus: the ARGO protocol. *Int Ophthalmol* 2023; 43(7):2315-2328
20. Roszkowska AM et al. Assessment of the predictive ability of theranostics for corneal cross-linking in treating keratoconus: a randomized clinical trial. *Ophthalmology* 2024;S0161-6420(24)00367-1

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**RegenSight**  
REGENERATE. IMPROVE

# *RegenSight*

*Delivering life-changing therapies to patients with keratoconus*

---